

eDNA isolation by precipitation.

Vs 4.2

17-08-2020

Materials:

-3M Sodium Acetate (NaAce) pH 5,2 -isopropanol -Phosphate Buffered Saline (PBS) -Purelink gDNA isolation kit (3957201, Brunschwig chemie bv) (190764-1L, Sigma-Aldrich) (15414499, Fischer Scientific) (10053293, Fischer Scientific)

- 1. Add 0.1 volume of 3 M Sodium Acetate to 1 volume of water sample (in a 15 or 50 ml tube).
- 2. Add 1 volume of isopropanol to the water sample.
- 3. Store overnight in 0-4 °C.
- 4. Centrifuge 1 h at 4600 x g in a centrifuge with swing-out rotor*.
- 5. Decant and discard the supernatant.
- 6. Centrifuge 5 min at 4600 x g.
- 7. Carefully pipet the supernatant and discard it.
- 8. Air-dry the pellet for 30 min at room temperature.
- 9. Resuspend the pellet in 200 µL 1x PBS.
- 10. Short spin down
- 11. Transfer the 200 μl to a 1,5 ml Eppendorf tube.
- 12. Store at 20°C or use directly in the Purelink gDNA isolation kit (step 13)

Use the following protocol to prepare lysate from the dissolved pellet. (Protocol extracted from manual Purelink gDNA isolation kit)

- 13. Set a water bath or heat block at 55°C.
- 14. Add 20 μ L Proteinase K (supplied with the kit) to the sample.
- 15. Add 20 μL RNase A (supplied with the kit) to the sample, mix well by brief vortexing and incubate at room temperature for 2 minutes.
- 16. Add 200 μL PureLink Genomic Lysis/Binding Buffer and mix well by vortexing to obtain a homogenous solution.
- 17. Incubate at 55°C for 10 minutes to promote protein digestion.
- 18. Spin for 5 min at 13000 x g and transfer the supernatant to a clean 1,5 ml Eppendorf tube.
- 19. Add 200 μL 96–100% ethanol to the lysate. Mix well by vortexing for 5 seconds to yield a homogenous solution.
- 20. Add the lysate (~640 μL) prepared with PureLink Genomic Lysis/Binding Buffer and ethanol to the PureLinkR Spin Column.
- 21. Centrifuge the column at 10,000 x g for 30 sec. at room temperature.
- 22. Discard the collection tube and place the spin column into a clean PureLink Collection Tube supplied with the kit.

- 23. Add 500 μL Wash Buffer 1 prepared with ethanol to the column.
- 24. Centrifuge column at room temperature at $10,000 \times g$ for 30 sec.
- 25. Discard the flow-through and place the spin column back into the PureLink collection tube.
- 26. Add 500 μ L Wash Buffer 2 prepared with ethanol to the column.
- 27. Centrifuge the column at 13000 x g for 2 minutes at room temperature. Discard collection tube.
- 28. Place the spin column into a new and clean PureLink collection tube.
- 29. Centrifuge the column at 13000 x g for 5 minutes at room temperature. Discard collection tube.
- 30. Place the PureLink Spin Column in an Eppendorf lowbind 1,5 ml ep.
- 31. Elute the DNA by adding 25 μ l elution buffer (provided in kit) and incubate for 1 minute at room temperature.
- 32. Centrifuge for 1 minute at 13000 x g.
- 33. The \approx 23 µl eluate, containing the eDNA, is ready to use in an qPCR assay.
- 34. For qPCR of point samples use duplicates of 5 µL purified sample

*do not spin above 4600 x g. The 15- or 50-ml tubes will break, although they are rated for 20 k x g this is not the case for a swing-out rotor where all the force concentrates in the tip of the tube!

Fixed angle rotor vs swingout rotor:

